



INTEGRATED MACHINE QA

myQA[®] Machines

Software, Detectors, Phantoms & Accessories

PROTECT +
ENHANCE +
SAVE LIVES



Machine QA ...

... Fundamental for Patient Safety



1 **Dosimetry QA**
Linac dose and system checked

2 **Imaging QA**
Linac imaging systems checked

3 **Morning QA**
Treatment system ready for the day

4 **Treatment Safety**
Peace of mind

Correct and accurate Linac Machine QA is a fundamental requirement for consistently safe and efficient patient treatments.

IBA Dosimetry is your proven partner to check and track all Machine QA needs, from daily to annual QA. Our imaging and dosimetry solutions provide highly efficient workflows and accurate analysis of your data.

IBA Dosimetry We Protect, Enhance and Save Lives.

At IBA we are passionate about providing innovative solutions for the diagnosis and treatment of cancer.

We focus on the well-being of patients, as well as the safe and efficient work of healthcare professionals worldwide.

IBA Dosimetry offers a full range of solutions for **Integrated Quality Assurance (QA)**, calibration procedures, and imaging markers, as well as services and training.

All our activities share a common goal: to **maximize efficiency and patient safety** in Radiation Therapy and Medical Imaging.



Innovation leadership in Machine QA

IBA Dosimetry's unique competence and leading innovations in Machine QA:

- ✓ > 1,500 satisfied customers worldwide trust IBA Dosimetry integrated Quality Assurance with myQA
- ✓ First complete Machine QA platform solution
- ✓ > 45 years of experience in Dosimetry and QA
- ✓ First Morning QA that combines efficiency with accuracy: myQA Daily
- ✓ 8 releases and software enhancements since the launch of myQA in 2015
- ✓ 24/7 support access from service teams in 3 time zones





myQA Machines, the complete protocol-based machine QA integrated on one platform.

- Full coverage of tests related to dosimetry, imaging, MLC QA... and more!
- Designed to integrate seamlessly with the myQA Platform
- Protocol-based machine QA (including TG-142 and other customizable protocols)
- Flexible scheduling tool to manage your tasks, resources, and time
- Integrate any of your hospital-specific tests with the individual tests module
- Interface to myQA Cockpit for quick and easy access to all QA results and trends
- Export any QA test result to a comprehensive report, traceable anytime
- Comprehensive analysis, archiving, and reporting tools

Integrated Software Platform



myQA[®] – The leading Integrated Quality Assurance

Over 1,500 healthcare centers in the world trust myQA as their global QA Platform.

myQA integrates myQA Machines to enable:

- End-to-end QA from commissioning to patient QA
- Machine QA data integration and comparison with beam scanning and patient QA
- All QA data on a central database in your hospital

myQA[®] Cockpit

Never Miss Any Key Information

- Have full confidence and control over all your Linacs, imaging devices and other technology
- In one view: verify which Machine QA tests have been completed or are still due for completion per schedule
- Detect negative trends before major problems arise with Trend Analysis



Dosimetry Plug-In

Perform automated dosimetry tests with the StarTrack, MatriXX, or Dolphin detectors, or by loading water phantom measurements.

Efficiency

- Acquire all key beam parameters in just one shot (dose output, profile analysis, energy verification)
- Analysis of main axes and diagonals (field size, symmetry, flatness, center, penumbra, light field)

1 Test Run ribbon for managing your tasks and their status. Finished tasks can be tracked in myQA Cockpit.

2 Detailed view of the reference and actual measurements (profiles).

3 Select your machine and visualize tasks due or add unscheduled tasks e.g. in case of maintenance.

4 Run your dosimetry test: Connect the detector and compare the actual measurement with your reference. Passes and fails are automatically displayed and the test status is recorded accordingly.

Direction	Parameter	Expected	Warning tolerance	Error tolerance	Actual	State
In-line	Center	-0.3 mm	± 1.5 mm	± 2.0 mm	-0.3 mm	Passed
In-line	Deviation	3.52 %	± 0.70 %	± 1.00 %	3.52 %	Passed
In-line	Flatness	108.72 %	± 1.00 %	± 1.00 %	108.72 %	Passed
In-line	Left Penumbra	9.3 mm	± 1.5 mm	± 2.0 mm	9.3 mm	Passed
In-line	Right Penumbra	8.9 mm	± 1.5 mm	± 2.0 mm	8.9 mm	Passed
In-line	Symmetry	100.85 %	± 0.70 %	± 1.00 %	100.85 %	Passed
Crossline	Center	1.0 mm	± 1.5 mm	± 2.0 mm	1.0 mm	Passed
Crossline	Deviation	2.83 %	± 0.70 %	± 1.00 %	2.83 %	Passed
Crossline	Flatness	102.86 %	± 0.70 %	± 1.00 %	102.86 %	Passed
Crossline	Left Penumbra	9.6 mm	± 1.5 mm	± 2.0 mm	9.6 mm	Passed
Crossline	Right Penumbra	9.6 mm	± 1.5 mm	± 2.0 mm	9.6 mm	Passed
Crossline	Symmetry	100.54 %	± 0.70 %	± 1.00 %	100.54 %	Passed

myQA FastTrack

Instant Measurements Plug-In

Real-time measurements & analysis with the StarTrack detector.

- Easy measurements outside scheduled routine checks
- Instant display of results, profile comparison and analysis (e.g. for beam steering, start-up behavior)
- Allows the Linac technician to set up unscheduled tests
- Import and export of measurements (ASCII)



MLC QA Plug-In

Automated MLC stripe tests ('picket fence test')

- EPID image analysis to determine leaf position accuracy & MLC transmission characteristics
- Identify if any MLC leaf is out of tolerance or which leaf number failed
- Verify your MLC at the 4 gantry cardinal angle
- Copy & paste detailed results for additional analysis

1 Intuitive and efficient: All MLC test parameters in one single interface view.

2 Test results at a glance: Review instantly if your MLC test passed, or see which test parameters exceeded your thresholds.

Baseline values	Expected	Warn +/-	Fail +/-	Actual
Line distance (mm)	50.14	0.20	0.50	50.14
Failing peaks (%)	0.00	0.30	0.50	0.00
Maximum difference (mm)	0.20	0.30	0.50	0.20
Standard deviation (mm)	0.07	0.30	0.50	0.07
Interstrip ratio (%)	113.01	3.00	5.00	113.01

3 Compare reference and actual MLC test. Overview of deviations. Adjust brightness and contrast easily.

4 Leaves deviation histogram: Easy and intuitive identification of the single stripe through color-coding.

VMAT QA Plug-In¹

Automated dynamic MLC QA and VMAT QA

- Automatic test analysis to verify accurate dose delivery using different dose rates, gantry speeds, and MLC leaf speeds
- Ensures that the changing rates and speeds during delivery do not adversely affect the delivered dose

1 Detailed numerical overview of the ROIs calculated characteristics.

2 Test results at a glance: Instant analysis of the correlation (ratio) between the open field and the specific delivered VMAT sequence.

1 ROI 1 mean: 0.996
ROI 2 mean: 1.003
ROI 3 mean: 1.003
ROI 4 mean: 0.997

2 ROI 1 standard deviation: 0.003
ROI 2 standard deviation: 0.001
ROI 3 standard deviation: 0.001
ROI 4 standard deviation: 0.003

ROI 1 name: 0.46 cm/s
ROI 2 name: 0.62 cm/s
ROI 3 name: 1.64 cm/s
ROI 4 name: 2.76 cm/s

New release

The New MLC Plug-In, for instant results and overview

- Complete Overview: All data on one screen incl. reference image and pass / fail display
- Completely automatic image pre-processing

2D-Imaging QA Plug-In

Automated imaging QA for EPIDs and for planar imaging [kV and MV]

- Fully automatic with all calculations performed in 5 seconds or less
- Compatible with all common imaging phantoms


New release

The New 2D-Imaging Plug-In, driving workflow efficiency

- Test efficiency: Pass / Fail clearly displayed
- Analysis power with the new special resolution MTF chart

1 Full overview: All EPID test parameters in one single interface view.

2 Test results at a glance: Instant overview of the CBCT test results and failed tests.



3 Reference and actual EPID test images. All regions of interest (RoI) are automatically detected.

4 Quick overview of the EPID's spatial resolution with the MTF display (Modulation Transfer Function).

3D-Imaging QA Plug-In

Automated imaging QA for CT and CBCT

- Including contrast, contrast to noise ratio, uniformity, HU deviation, spatial resolution, imaging scaling, and more!
- Compatible with all common imaging phantoms


New release

The New 3D-Imaging Plug-In, now enhanced and even faster

- Save time: Complete overview on one screen
- Instant Confirmation: Clear display of pass / fail
- Workflow ease of use: From image import to the new MTF chart

1 Intuitive & efficient: All CBCT test parameters in one single interface view.

2 Test results at a glance: Instant overview of the CBCT test results and failed tests.



3 Reference and actual CT / CBCT test image. Automatic detection and analysis of Regions of Interest.

4 Instant overview of the spatial resolution of your CBCT system with the MTF display (Modulation Transfer Function).



Integrate Any Of Your Specific Tests In myQA Machines

Your individual QA checks and any other tests are easily integrated with the myQA customizable generic tests functionality


- Easily compose any test you need to check and document
- Import of existing data and tests from Excel
- Make myQA your complete solution for:
 - everything you need to check regularly
 - all tests you need to track and record in one database
- Examples
 - check the temperature of your medicine fridge
 - check your room lighting or security locks ...

With the myQA's Individual Test feature we even integrate and track checks such as "Doors Locked" or "Oxygen Off".



James P Nunn, MS, CHP, DABR
Senior Medical Physicist,
LewisGale Regional
Cancer Centre,
Pulaski, USA

Detectors for Machine QA



myQA Daily – Fast and easy wireless setup and morning QA test results

Smartly designed measurement tools are your basis for efficient and accurate Linac QA.

IBA Dosimetry offers a wide range of dedicated solutions to make your daily, monthly, and annual QA the fastest, most accurate, and most reliable.


- Designed to integrate seamlessly with myQA
- Robust for long lasting performance
- Accuracy based on ionization chamber technology

New release

myQA Daily

The only solution for Fast, Easy, and Accurate Morning Linac QA.

- The largest number of ionization chambers (125) of any morning QA detector provides more beam data for more accurate beam quality verification
- The web browser-based myQA Daily application allows flexible test execution from any network PC or tablet and easy access to test results
- Compatible with Varian Halcyon™ and all standard linacs
- Integrated energy checks



The software provides an intuitive overview of the accuracy of all tests. Each test result can be verified in a detailed view compared with the expected result.

StarTrack

Your High-End Detector for Advanced Machine QA


- All main tests in one shot: dose, profiles, diagonals energy verification, etc.
- 453 air-vented ionization chambers with optimized geometry for Machine QA
- Convenient beam constancy verification in one single shot using specific build-up plates
- Automatic k(t,p) correction
- Parallel readout from independent electrometers
- Instant results and real-time analysis using the Dosimetry plug-in for myQA Machines
- Patented energy verification method
- Tabletop or gantry mount (optional)



MatriXX

Your Flexible Detector


- Detector for Patient QA (with myQA Patients software) as well as for Machine QA
- Connect your MatriXX to myQA Machines for fast and accurate Linac Machine QA
- Choose the right detector from the MatriXX family (MatriXX Evolution / MatriXX FFF)
- Patented energy verification method
- Tabletop or gantry mount (optional)



Build-Up Plates

For Energy Constancy Verification

- Specific build-up plates for StarTrack and MatriXX detectors
- Convenient beam constancy verification in one single shot




Gantry Mount

- To detect dosimetry errors introduced with Linac rotation angles
- Available for MatriXX and StarTrack to attach your detector to all major linac accessory mount interfaces



Phantoms for Machine QA



Integrated phantoms and test devices for your comprehensive Machine and Imaging QA needs.

Select from a range of phantoms for dedicated machine QA tasks.

- IBA Dosimetry imaging phantoms for 2D and 3D image quality verification
- Integrate your existing phantoms (myQA supports most common imaging phantoms)




myQA has given me full control of my data by connecting all QA applications on one platform and into one central database. With myQA, the quality assurance becomes schedulable - in every sense of the word. Another highlight for me is the web-based myQA Cockpit dashboard which allows us to quickly retrieve our machine QA status updates anywhere in the department. myQA is truly an all-in-one solution.

Luis Brualla González
Hospital General Universitario,
ERESA, Valencia, Spain

X-LITE

Light Field Alignment Checks

- Fast and precise check of the radiation and light fields without film or additional hardware
- Easy setup against the light field
- Fluorescent plate visualizes your radiation field right after the irradiation without additional processing
- Field scales marked with 5 × 5, 10 × 10, 15 × 15, and 20 × 20 cm²



Cylindrical Phantom

Dose Constancy and Isocenter Check

- Verify the mechanical stability of gantry/imager position (CBCT/ EPID) with a small steel ball insert (Winston-Lutz Test)
- Measurement of dose constancy in various gantry angles and in rotational beams
- Adaptors available for most common ion chambers

For more information please refer to the IBA white paper 'Tg-142_Daily Generic Tests' and 'Winston-Lutz & Star Shot Test'




Example of a verification image using the IBA Cylindrical Phantom.

Disk Phantom

Isocenter Verification with Film

- Easy and precise method of verifying isocenter accuracy (e.g. for stereotactic applications, star shot)
- Isocenter is determined by an appendant marker




IBA Dosimetry Disk Phantom: Film with a typical star shot pattern.

Additional QA Hardware Available

- Full set of chambers and diodes
- Reference class electrometers
- Plastic slab Phantoms and chamber inserts
- Round CT and RTPS Phantom



For more information visit iba-dosimetry.com



Imaging QA in Radiation Therapy



1 Setup test plate & MagicMaX XR Multi-Detector for QA of kV imager



2 Comprehensive instant analysis of imaging dose

3 Reporting in myQA

Image Quality & Imaging Dose QA


The affordable solution for all your x-ray and CT imaging QA needs for image quality and dose. Also available as complete solution kits.

- TG-142: Kit provides all the tools necessary for the Imaging Dose verification with Multimeter MagicMaX
- CyberKnife Imaging QA: Unique efficiency through support of serial exposures

Primus A

Test plate for kV planar image QA

- Easy image QA of your IGRT imaging systems or flat-panel imager (EPID)
- Verify complete contrast determination, special resolution, scaling discrepancy, uniformity and positioning offsets



CT Phantom

- Innovative 3-part nested PMMA phantom for CTDI measurements
- Designed to image pediatric and adult head and body
- According to FDA performance standard for diagnostic X-ray systems (21CFR 1020.33)




MagicMaX

Imaging-Dose Multimeter

- Fast, simple, and accurate beam analysis and dosimetry for your CBCT, OBI, and CT-Sim, and 2D/3D imaging
- Ideally suited for Varian OBI, Elekta or CyberKnife
- In a single exposure, evaluate your kV beam and imaging dose or flat-panel imager (with the Primus L Test Plate)
- MagicMaX with exchangeable detectors for kV or CT dose measurements



The Flexible Solution For All Imaging QA Needs In The RT Department



CT virtual Sim

Varian OBI kV /CBCT




Elekta kV /CBCT

CyberKnife kV

Technical Specifications

StarTrack Detector	
Energy Range	Photons: ^{60}Co , 4–18 MV, flattened and FFF beams. Electrons: 6–21 MeV.
Dose Linearity	0.5 % from 10 cGy to 5 Gy integral dose. 0.5 % from 0.1 Gy/min up to 4 Gy/min dose rate.
k_{Tp} Correction	Temperature [10–40 °C], pressure [70–110 kPa].
Sensor Layout	Chamber arrays organized along main axes and diagonals, 8 additional chambers for energy constancy check.
Spatial Resolution	5mm for horizontal and vertical lines. 7mm for diagonals.
Chamber Type	Vented plane parallel ionization chambers.
Chamber Size	Cylindrical, 3 [Ø] x 5 [h] mm, sensitive volume 35 mm ³ .
Typical Sensitivity	1.1 nC/Gy [^{60}Co]
Electrometer	8 TERA ASICs [each contains 64 independent electrometers].
Sampling Time	min. 10 ms
Readout	Parallel and synchronous readout with no dead time.
myQA Daily Detector	
Energy Range	Photons: ^{60}Co , 4–24 MV, flattened and FFF beams. Electrons: 4–24 MeV
Dose Linearity	<0,5 % for dose greater than 0,25 Gy and dose rate greater than 0,3 Gy/min
k_{Tp} Correction	Yes
Sensor Layout	125 ionization chambers, layout optimized for 10×10 and 20×20 cm ² field measurements
Spatial Resolution	5 mm grid
Chamber Type	Vented plane parallel ionization chambers
Chamber Size	3,2 mm Ø, 2 mm height, volume 16 mm ³
Typical Sensitivity	0,53 nC/Gy [^{60}Co]
Electrometer	Integrated 128 channel electrometer chip
Sampling Time	500 ms
Readout	Parallel and synchronous readout with no dead time.

Follow us

-  YouTube | youtube.com/user/ibadosimetry
-  LinkedIn | linkedin.com/company/iba-dosimetry-gmbh
-  Twitter | twitter.com/ibadosimetry

Machine-QA-Br-E_Rev1_0719 | © IBA 2019 | All rights reserved
Depicted product images may differ from the actual scope of delivery.
Images and technical specifications are subject to change without prior notice.

IBA Dosimetry

Integrated Quality Assurance
Europe, Middle East, Africa | +49-9128-6070
North America and Latin America | +1-901-386-2242
Asia Pacific | +86-10-8080-9288
dosimetry-info@iba-group.com | iba-dosimetry.com

